

Nuclear factor NF-kappa-B p105 subunit(NFKB1)

OVERVIEW

PRODUCT REFERENCE	PX-P4764
ORIGIN SPECIES	
FRAGMENT	Asp40-Tyr349
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	35.98kDa
FORM	
BUFFER	PBS PH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P19838
ALIAS/SYNONYMS	DNA-binding factor KBF1,EBP-1,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1,