

Lithostathine-1-beta(REG1B)

OVERVIEW

PRODUCT REFERENCE	PX-P4781
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln23-Asn166
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.23 kDa
FORM	
BUFFER	PBS pH 7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P48304
ALIAS/SYNONYMS	Pancreatic stone protein 2,PSP-2,Regenerating islet-derived protein 1-beta,REG-1-beta,Regenerating protein I beta,PSPS2, REGL

