

Interleukin-31 receptor subunit alpha(IL31RA)

OVERVIEW

PRODUCT REFERENCE	PX-P4802
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Val130-Asp622
TAG	N terminus His tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	55.75 kDa
FORM	
BUFFER	PBS pH 7.5, 0.02% SKL
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8NI17
ALIAS/SYNONYMS	IL-31 receptor subunit alpha, IL-31R subunit alpha, IL-31R-alpha, IL-31RA, Cytokine receptor-like 3, GLM-R, hGLM-R, Gp130-like monocyte receptor, Gp130-like receptor, ZcytoR17, CRL3, GPL

