

Meiotic recombination protein REC8 homolog(Rec8)

OVERVIEW

PRODUCT REFERENCE	PX-P4828
ORIGIN SPECIES	
FRAGMENT	Met1-Pro591
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	67.42kDa
FORM	
BUFFER	TBS, pH8.0
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8C5S7
ALIAS/SYNONYMS	Cohesin Rec8p, Mei8, Rec8L1

