

Nuclear receptor ROR-alpha(RORA)

OVERVIEW

PRODUCT REFERENCE	PX-P4834
ORIGIN SPECIES	
FRAGMENT	Glu263-Gly523
TAG	N-terminal His Tag & sumo Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	42,9kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P35398
ALIAS/SYNONYMS	Nuclear receptor RZR-alpha, Nuclear receptor subfamily 1 group F member 1, RAR-related orphan receptor A, Retinoid-related orphan receptor-alpha, NR1F1, RZRA

