

Isoform SREBP-1C of Sterol regulatory element-binding protein 1

OVERVIEW

PRODUCT REFERENCE	PX-P4844
ORIGIN SPECIES	
FRAGMENT	Met1-Leu466
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	50.2kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P36956-3
ALIAS/SYNONYMS	SREBP-1, Class D basic helix-loop-helix protein 1, bHLHd1, Sterol regulatory element-binding transcription factor 1, BHLHD1, SREBP1