

Indoleamine 2,3-dioxygenase 1(IDO1)

OVERVIEW

PRODUCT REFERENCE	PX-P4856
ORIGIN SPECIES	
FRAGMENT	Met1-Gly403
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	47.6kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P14902
ALIAS/SYNONYMS	IDO-1, Indoleamine-pyrrole 2,3-dioxygenase, IDO, INDO

