

Interleukin-5(IL5)

OVERVIEW

PRODUCT REFERENCE	PX-P4859
ORIGIN SPECIES	
FRAGMENT	Ile20-Ser134
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.3kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05113
ALIAS/SYNONYMS	IL-5,B-cell differentiation factor I,Eosinophil differentiation factor,T-cell replacing factor,TRF

