

Lysostaphin(Iss)

OVERVIEW

PRODUCT REFERENCE	PX-P4876
ORIGIN SPECIES	
FRAGMENT	Ala248-Lys493
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	35.97kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P10547
ALIAS/SYNONYMS	Glycyl-glycine endopeptidase

