

Inhibin alpha chain

OVERVIEW

PRODUCT REFERENCE	PX-P4882
ORIGIN SPECIES	
FRAGMENT	Ser233-Ile366
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.35kDa
FORM	
BUFFER	
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05111
ALIAS/SYNONYMS	A-inhibin subunit; INHA; Inhibin A; inhibin alpha chain; inhibin, alpha

